

TRIUMPH TIGER 800 /XC/XCA/XCX/XR/XRT/XRX '11-'19
TRIUMPH TIGER 900 /GT '20

M O U N T I N G I N S T R U C T I O N S

REF 20383N

PROTECTOR BASCULANTE PHB19 / PHB19 PROTECTOR SWINGARM



1 Person



15 min



EASY



TUBE





Nº 10

RATCHET



Nº 10

ID.	PIEZA / PART	QTY.	REF.:
1		1	PROBAS33-3-0033
2		1	PROBAS22-3-0022
3		1	VARIFO14-0-0014
4		2	TUZBM6000000024
5		2	ARANDEPD125-054
6		2	TAPAFINAL17-18


TRIUMPH TIGER 800 /XC/XCA/XCX/XR/XRT/XRX '11-'19 TRIUMPH TIGER 900 /GT '900

REF 20383N

INSTRUCCIONES DE MONTAJE - MOUNTING INSTRUCTIONS - MONTAGEANLEITUNG


1

ASSEMBLY
MONTAJE




2

ASSEMBLY
MONTAJE




3

ASSEMBLY
MONTAJE



CONSULTAR EN NUESTRA PAGINA WEB LA
DISPONIBILIDAD DEL VIDEO DE MONTAJE DE
ESTA REFERENCIA.
CONSULT OUR WEBSITE THE VIDEO MOUNT
AVAILABILITY OF THIS REFERENCE.



6 Extracción del tapón / Removing the cap